

Nivatrotamab

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Nivatrotamab (Hu3F8-BsAb), a humanized anti-GD2/CD3 bispecific antibody and T-cell engager specific to CD3 and GD2, is employed in the study of neuroblastoma [1].
Targets(IC50)	Immunology/Inflammation related,Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481