

Obrindatamab

Chemical Properties

CAS No. : 2069959-72-2

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Obrindatamab, a humanized anti-B7-H3/CD3 bispecific antibody, targets B7-H3 and CD3, facilitating redirected cytotoxic T-lymphocyte (CTL) activity towards B7-H3-expressing cancer cells. It is utilized in cancer research [1].
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481