

Ongericimab

Chemical Properties

CAS No. : 2145096-91-7

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ongericimab (JS002) is a recombinant humanized monoclonal antibody targeting proprotein convertase subtilisin/kexin type 9 (PCSK9). By blocking the binding of PCSK9 to the low-density lipoprotein receptor (LDLR), the antibody reduces LDLR degradation, thereby increasing hepatic LDL-C uptake and lowering plasma cholesterol levels. It is primarily used to research hypercholesterolemia and cardiovascular risk management.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481