

## Praluzatamab

## Chemical Properties

CAS No. :	2145109-70-0
Formula:	
Molecular Weight:	146.9 kDa
Storage:	Store at low temperature -20°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Praluzatamab is a monoclonal antibody targeting an activated white blood cell adhesion molecule (ALCAM/CD116). It is a naked antibody used to synthesize antibody active molecule conjugates (Praluzatamab ravnansine) and can be employed to study breast cancer.
Targets(IC50)	ADC Antibody, Immunology/Inflammation related

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0068 mL	0.034 mL	0.0681 mL
5 mM	0.0014 mL	0.0068 mL	0.0136 mL
10 mM	0.0007 mL	0.0034 mL	0.0068 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481