

Pacanalotamab

Chemical Properties

CAS No. : 2251756-52-0

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pacanalotamab (AMG 420; BI-836909) is a bispecific T-cell engager (BiTE) that specifically targets BCMA (B cell maturation antigen) and CD3. By redirecting T cells to cells expressing BCMA, Pacanalotamab facilitates the T-cell mediated lysis of human multiple myeloma (MM) cell lines [1].
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481