

Quisovalimab

Chemical Properties

CAS No. : 2427667-03-4

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Quisovalimab is a fully human monoclonal antibody that specifically binds to and neutralizes LIGHT(TNFSF14), a member of tumor necrosis factor superfamily, thus blocking its interaction with the receptor and inhibiting T cell-related inflammatory response. Quisovalimab has been used to study the inflammatory pathway mediated by LIGHT, and it has potential therapeutic value in COVID-19-related acute respiratory distress syndrome (ARDS) and cytokine release syndrome (CRS).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481