

Pavurutamab

Chemical Properties

CAS No. : 2250292-39-6

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pavurutamab (AMG-701) is a bispecific T-cell engager targeting CD3 and B-cell maturation antigens (BCMA), with an extended half-life derived from Pacanalotamab. Its Fc region is engineered for improved pharmacokinetics, offering potential in immune regulation and treatment of multiple myeloma (MM) [1] [2] [3] [4].
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481