

Pegaldesleukin

Chemical Properties

CAS No. : 141752-91-2

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pegaldesleukin, a conjugate of polyethylene glycol and interleukin-2 (PEG-IL2), exhibits antiviral activity with potential applications in managing HIV. It may delay the progression of HIV infection by preserving the immune repertoire [1] [2] [3].
Targets(IC50)	HIV Protease,Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481