

## Regavirumab

## Chemical Properties

CAS No. : 153101-26-9

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Regavirumab (MCA C23), a human monoclonal antibody, targets cytomegalovirus (CMV), human serum albumin (HSA), and amino acetic acid [1].
Targets(IC50)	Virus Protease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481