

Dalutrafusp alfa

Chemical Properties

CAS No. : 2419918-89-9

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dalutrafusp alfa (AGEN-1423; GS-1423) is a bifunctional antibody targeting CD73 and TGF- β , components of the immunosuppressive pathway [1].
Targets(IC50)	CD73,TGF-beta/Smad

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481