

Pexelizumab

Chemical Properties

CAS No. : 219685-93-5

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pexelizumab (h5G1.1-SC), a humanized scFv monoclonal antibody, targets the C5 complement component, thereby inhibiting apoptosis and leukocyte infiltration. It is utilized in researching cerebral ischemia-reperfusion (IR) injury and myocardial infarction [1] [2].
Targets(IC50)	Apoptosis, Complement System

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481