

Pimivalimab

Chemical Properties

CAS No. : 2293951-22-9

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pimivalimab (JTX-4014) is a high-affinity, fully human IgG4 monoclonal antibody targeting the cell surface receptor PD-1. By blocking the interaction between PD-1 and its ligands PD-L1 and PD-L2, the antibody reverses the immunosuppressed state of T cells and enhances the anti-tumor immune response against various solid tumors. It is frequently studied as a monotherapy or in combination with other therapies, such as LILRB2 inhibitors.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481