

## Sevirumab

## Chemical Properties

CAS No. : 138660-96-5  
Formula:  
Molecular Weight:  
Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Sevirumab (MSL-109), a human IgG1 monoclonal antibody, neutralizes cytomegalovirus (CMV) by specifically targeting CMV gH complexes, effectively inhibiting viral replication with an EC50 value of 0.3 µg/mL [1] [2].
Targets(IC50)	Virus Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481