

## Efineptakin alfa

## Chemical Properties

CAS No. : 2026634-47-7

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Efineptakin alfa (NT-17), a long-acting recombinant human interleukin-7 (IL-7), promotes the proliferation and survival of CD4+ and CD8+ cells in both humans and mice. It has applications in glioblastoma research [1].
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481