

Eflapegrastim

Chemical Properties

CAS No. : 1384099-30-2

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Eflapegrastim, a humanized IgG4 monoclonal antibody and granulocyte colony-stimulating factor (G-CSF), targets the G-CSF receptor (c-Fms). It promotes the proliferation and differentiation of neutrophil progenitor cells, ensuring a stable count of mature and functional neutrophils, and effectively shortens the duration of neutropenia [1].
Targets(IC50)	c-Fms

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481