

Efruxifermin

Chemical Properties

CAS No. : 2375240-92-7

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Efruxifermin, an Fc-FGF21 fusion protein consisting of a human IgG1 Fc domain linked to modified human FGF21, exhibits a prolonged half-life and enhanced receptor affinity compared to its native counterpart. It is utilized in researching non-alcoholic steatohepatitis [1].
Targets(IC50)	FGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481