


Sialyl lewis x-lactose

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 127923-85-7 |
| Formula: | C43H72N2O33 |
| Molecular Weight: | 1145.03 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|---|
| Description | Sialyl Lewis X-Lactose, a biochemical reagent, serves as a biological material or organic compound for life science research. |
| Targets(IC50) | Others,Endogenous Metabolite |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 0.8733 mL | 4.3667 mL | 8.7334 mL |
| 5 mM | 0.1747 mL | 0.8733 mL | 1.7467 mL |
| 10 mM | 0.0873 mL | 0.4367 mL | 0.8733 mL |
| 50 mM | 0.0175 mL | 0.0873 mL | 0.1747 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481