


TBC-3486

Chemical Properties

CAS No. : 247044-77-5
 Formula: C₂₇H₃₁N₃O₆S₂
 Molecular Weight: 557.68
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	TBC-3486 is a novel potent inhibitor of integrin α 4, being 200-fold more potent in inhibiting α 4 β 1 than α 4 β 7.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7931 mL	8.9657 mL	17.9314 mL
5 mM	0.3586 mL	1.7931 mL	3.5863 mL
10 mM	0.1793 mL	0.8966 mL	1.7931 mL
50 mM	0.0359 mL	0.1793 mL	0.3586 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481