

## Compound T128317

## Chemical Properties


CAS No. :

Formula: C27H32O14

Molecular Weight: 580.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

2-(3,4-dimethoxyphenyl)-5,6,7,8-tetramethoxy-3-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)-4H-chromen-4-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7225 mL	8.6127 mL	17.2253 mL
5 mM	0.3445 mL	1.7225 mL	3.4451 mL
10 mM	0.1723 mL	0.8613 mL	1.7225 mL
50 mM	0.0345 mL	0.1723 mL	0.3445 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481