

## 3-Epimeliasenin B

## Chemical Properties


CAS No. : 1222475-77-5

Formula: C<sub>30</sub>H<sub>44</sub>O<sub>4</sub>

Molecular Weight: 468.678

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

3-Epimeliasenin B is a natural product for research related to life sciences. The catalog number is TN5953 and the CAS number is 1222475-77-5.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1337 mL	10.6683 mL	21.3365 mL
5 mM	0.4267 mL	2.1337 mL	4.2673 mL
10 mM	0.2134 mL	1.0668 mL	2.1337 mL
50 mM	0.0427 mL	0.2134 mL	0.4267 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481