

Benzenesulfonic acid, 3-(2-(4-(2-(4-((2-cyanoethyl)ethylamino)phenyl)diazenyl)-3-methoxyphenyl)diazenyl)-


Chemical Properties

CAS No. : 68227-58-7

Formula: C₂₄H₂₄N₆O₄S

Molecular Weight: 492.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzenesulfonic acid, 3-(2-(4-(2-(4-((2-cyanoethyl)ethylamino)phenyl)diazenyl)-3-methoxyphenyl)diazenyl)- is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0303 mL	10.1513 mL	20.3025 mL
5 mM	0.4061 mL	2.0303 mL	4.0605 mL
10 mM	0.203 mL	1.0151 mL	2.0303 mL
50 mM	0.0406 mL	0.203 mL	0.4061 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481