


Polygalasaponin XXXV

Chemical Properties

CAS No. :	184479-28-5
Formula:	C63H98O31
Molecular Weight:	1351.43
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Polygalasaponin XXXV is a natural product for research related to life sciences. The catalog number is TN6533 and the CAS number is 184479-28-5.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.740 mL	3.6998 mL	7.3996 mL
5 mM	0.148 mL	0.740 mL	1.4799 mL
10 mM	0.074 mL	0.370 mL	0.740 mL
50 mM	0.0148 mL	0.074 mL	0.148 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481