


3,4-Dihydro-6,7-(methylenedioxy)-2(1H)-quinolinone

Chemical Properties

CAS No. :	94527-34-1
Formula:	C ₁₀ H ₉ NO ₃
Molecular Weight:	191.186
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3,4-Dihydro-6,7-(methylenedioxy)-2(1H)-quinolinone is a useful organic compound for research related to life sciences. The catalog number is T123884 and the CAS number is 94527-34-1.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2304 mL	26.152 mL	52.304 mL
5 mM	1.0461 mL	5.2304 mL	10.4608 mL
10 mM	0.523 mL	2.6152 mL	5.2304 mL
50 mM	0.1046 mL	0.523 mL	1.0461 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481