

HAE acetate(64111-99-5 free base)

Chemical Properties

CAS No. :

Formula: C16H25N5O8

Molecular Weight: 415.4

His-Ala-Glu(acetate salt)

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HAE acetate is a 3-amino acid peptide composed of histidine, alanine and glutamate
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4073 mL	12.0366 mL	24.0732 mL
5 mM	0.4815 mL	2.4073 mL	4.8146 mL
10 mM	0.2407 mL	1.2037 mL	2.4073 mL
50 mM	0.0481 mL	0.2407 mL	0.4815 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481