


N.0096

## Chemical Properties

CAS No. : 75052-76-5  
 Formula: C<sub>28</sub>H<sub>23</sub>ClNaO<sub>4</sub>P  
 Molecular Weight: 512.9  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N.0096 is a thromboxane synthesis inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9497 mL	9.7485 mL	19.497 mL
5 mM	0.3899 mL	1.9497 mL	3.8994 mL
10 mM	0.195 mL	0.9748 mL	1.9497 mL
50 mM	0.039 mL	0.195 mL	0.3899 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481