


Olcorolimus

Chemical Properties

CAS No. :	186752-78-3
Formula:	C ₅₁ H ₈₁ N ₁ O ₁₂
Molecular Weight:	900.19
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Olcorolimus inhibits single allergen-induced allergic inflammation such as eosinophilic and lymphocytic infiltration and mRNA expression for interleukin (IL)-4 and IL-5
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1109 mL	5.5544 mL	11.1088 mL
5 mM	0.2222 mL	1.1109 mL	2.2218 mL
10 mM	0.1111 mL	0.5554 mL	1.1109 mL
50 mM	0.0222 mL	0.1111 mL	0.2222 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481