


LS 105790

Chemical Properties

CAS No. : 137734-07-7
 Formula: C72H95N3O19
 Molecular Weight: 1306.53
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	LS 105790 is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7654 mL	3.8269 mL	7.6539 mL
5 mM	0.1531 mL	0.7654 mL	1.5308 mL
10 mM	0.0765 mL	0.3827 mL	0.7654 mL
50 mM	0.0153 mL	0.0765 mL	0.1531 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481