

## Compound T128838

## Chemical Properties

CAS No. :

Formula: C<sub>41</sub>H<sub>68</sub>O<sub>11</sub>

Molecular Weight: 736.98

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	2-({1-[(4E)-2-hydroxy-6-methylhepta-4,6-dien-2-yl]-3a,3b,6,6,9a-pentamethyl-7-[(3,4,5-trihydroxyoxan-2-yl)oxy]-hexadecahydro-1H-cyclopenta[a]phenanthren-11-yl}oxy)-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3569 mL	6.7844 mL	13.5689 mL
5 mM	0.2714 mL	1.3569 mL	2.7138 mL
10 mM	0.1357 mL	0.6784 mL	1.3569 mL
50 mM	0.0271 mL	0.1357 mL	0.2714 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481