


5-Tricosyl-1,3-benzenediol

Chemical Properties

CAS No. :	70110-60-0
Formula:	C ₂₉ H ₅₂ O ₂
Molecular Weight:	432.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Tricosyl-1,3-benzenediol is a natural product for research related to life sciences. The catalog number is TN3138 and the CAS number is 70110-60-0.
Targets(IC50)	Others,DNA/RNA Synthesis

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.311 mL	11.5548 mL	23.1096 mL
5 mM	0.4622 mL	2.311 mL	4.6219 mL
10 mM	0.2311 mL	1.1555 mL	2.311 mL
50 mM	0.0462 mL	0.2311 mL	0.4622 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481