


## Maropitant citrate anhydrous

## Chemical Properties

CAS No. :	862543-54-2
Formula:	C <sub>38</sub> H <sub>48</sub> N <sub>2</sub> O <sub>8</sub>
Molecular Weight:	660.8
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Maropitant citrate anhydrous is a neurokinin (NK1) receptor antagonist.
Targets(IC50)	Others, Neurokinin receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5133 mL	7.5666 mL	15.1332 mL
5 mM	0.3027 mL	1.5133 mL	3.0266 mL
10 mM	0.1513 mL	0.7567 mL	1.5133 mL
50 mM	0.0303 mL	0.1513 mL	0.3027 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481