

2-hydroxy-3-[5-(methoxycarbonyl)cyclohexa-2,4-dien-1-yl]propanoic acid

Chemical Properties


CAS No. :

Formula: C₁₁H₁₄O₅

Molecular Weight: 226.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-hydroxy-3-[5-(methoxycarbonyl)cyclohexa-2,4-dien-1-yl]propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4203 mL	22.1014 mL	44.2028 mL
5 mM	0.8841 mL	4.4203 mL	8.8406 mL
10 mM	0.442 mL	2.2101 mL	4.4203 mL
50 mM	0.0884 mL	0.442 mL	0.8841 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481