

Taikuguasin D aglycon

Chemical Properties


CAS No. : 2185847-09-8

Formula: C₃₁H₅₀O₄

Molecular Weight: 486.737

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

Taikuguasin D aglycon is a natural product for research related to life sciences. The catalog number is TN6311 and the CAS number is 2185847-09-8.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0545 mL	10.2724 mL	20.5448 mL
5 mM	0.4109 mL	2.0545 mL	4.109 mL
10 mM	0.2054 mL	1.0272 mL	2.0545 mL
50 mM	0.0411 mL	0.2054 mL	0.4109 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481