

Daphnodorin B

Chemical Properties


CAS No. :

Formula: C₃₀H₂₂O₁₀

Molecular Weight: 542.496

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Daphnodorin B is a useful organic compound for research related to life sciences and the catalog number is T123854.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8433 mL	9.2166 mL	18.4332 mL
5 mM	0.3687 mL	1.8433 mL	3.6866 mL
10 mM	0.1843 mL	0.9217 mL	1.8433 mL
50 mM	0.0369 mL	0.1843 mL	0.3687 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481