


2'',3''-Dihydro-3',3'''-biapigenin

Chemical Properties

CAS No. :	151455-25-3
Formula:	C ₃₀ H ₂₀ O ₁₀
Molecular Weight:	540.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'',3''-Dihydro-3',3'''-biapigenin is a natural product for research related to life sciences. The catalog number is TN5496 and the CAS number is 151455-25-3.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8502 mL	9.2512 mL	18.5024 mL
5 mM	0.370 mL	1.8502 mL	3.7005 mL
10 mM	0.185 mL	0.9251 mL	1.8502 mL
50 mM	0.037 mL	0.185 mL	0.370 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481