


Tolmetin β -D-Glucuronide

Chemical Properties

CAS No. :	71595-19-2
Formula:	C ₂₁ H ₂₃ N _O ₉
Molecular Weight:	433.41
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tolmetin β -D-Glucuronide has reactive properties with human serum albumin, allowing identification of the binding site and reaction mechanism by tandem mass spectrometry.
Targets(IC50)	Others, Drug Metabolite

Solubility Information

Solubility	DMSO: Slightly soluble Methanol: Slightly soluble, Heating is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3073 mL	11.5364 mL	23.0728 mL
5 mM	0.4615 mL	2.3073 mL	4.6146 mL
10 mM	0.2307 mL	1.1536 mL	2.3073 mL
50 mM	0.0461 mL	0.2307 mL	0.4615 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481