

5-[(2Z,8Z,11Z)-pentadeca-2,8,11-trien-1-yl]benzene-1,3-diol

Chemical Properties


CAS No. :

Formula: C₂₁H₃₀O₂

Molecular Weight: 314.47

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-[(2Z,8Z,11Z)-pentadeca-2,8,11-trien-1-yl]benzene-1,3-diol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.180 mL	15.8998 mL	31.7995 mL
5 mM	0.636 mL	3.180 mL	6.3599 mL
10 mM	0.318 mL	1.590 mL	3.180 mL
50 mM	0.0636 mL	0.318 mL	0.636 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481