

Bis-PEG15-acid

Chemical Properties

CAS No. :

Formula: C34H66O19

Molecular Weight: 778.88

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bis-PEG15-acid, a PEG-based linker for PROTACs, joins two essential ligands crucial for forming PROTAC molecules and enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binds to an E3 ubiquitin ligase, and the other to the target protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2839 mL	6.4195 mL	12.8389 mL
5 mM	0.2568 mL	1.2839 mL	2.5678 mL
10 mM	0.1284 mL	0.6419 mL	1.2839 mL
50 mM	0.0257 mL	0.1284 mL	0.2568 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481