

N-Me-N-bis(PEG4-acid)

Chemical Properties

CAS No. : 2055014-77-0


Formula: C23H45NO12

Molecular Weight: 527.6

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N-Me-N-bis(PEG4-acid) is a PEG-based linker used in PROTACs to join two essential ligands, facilitating the formation of PROTAC molecules. This linker enables selective protein degradation by utilizing the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets the desired protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8954 mL	9.4769 mL	18.9538 mL
5 mM	0.3791 mL	1.8954 mL	3.7908 mL
10 mM	0.1895 mL	0.9477 mL	1.8954 mL
50 mM	0.0379 mL	0.1895 mL	0.3791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481