

Compound N085-0035

Chemical Properties


CAS No. :

Formula: C₁₃H₂₄NO₂

Molecular Weight: 226.339

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N085-0035 is a useful organic compound for research related to life sciences and the catalog number is T131168.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4181 mL	22.0907 mL	44.1813 mL
5 mM	0.8836 mL	4.4181 mL	8.8363 mL
10 mM	0.4418 mL	2.2091 mL	4.4181 mL
50 mM	0.0884 mL	0.4418 mL	0.8836 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481