


## Saikochromone A

## Chemical Properties

CAS No. :	132624-99-8
Formula:	C <sub>11</sub> H <sub>10</sub> O <sub>5</sub>
Molecular Weight:	222.19
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Saikochromone A is a natural product for research related to life sciences. The catalog number is TN6316 and the CAS number is 132624-99-8.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5007 mL	22.5033 mL	45.0065 mL
5 mM	0.9001 mL	4.5007 mL	9.0013 mL
10 mM	0.4501 mL	2.2503 mL	4.5007 mL
50 mM	0.090 mL	0.4501 mL	0.9001 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481