


Kelthane-S

Chemical Properties

CAS No. : 52869-78-0
 Formula: C₁₉H₂₁Cl₅NO₅PS
 Molecular Weight: 583.68
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Kelthane-S is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7133 mL	8.5663 mL	17.1327 mL
5 mM	0.3427 mL	1.7133 mL	3.4265 mL
10 mM	0.1713 mL	0.8566 mL	1.7133 mL
50 mM	0.0343 mL	0.1713 mL	0.3427 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481