

## Compound N050-0025

## Chemical Properties


CAS No. :

Formula: C<sub>31</sub>H<sub>52</sub>O<sub>2</sub>

Molecular Weight: 456.755

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N050-0025 is a useful organic compound for research related to life sciences and the catalog number is T131244.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1893 mL	10.9467 mL	21.8933 mL
5 mM	0.4379 mL	2.1893 mL	4.3787 mL
10 mM	0.2189 mL	1.0947 mL	2.1893 mL
50 mM	0.0438 mL	0.2189 mL	0.4379 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481