

## Palbociclib-Succinic acid

## Chemical Properties


CAS No. :

Formula: C28H33N7O5

Molecular Weight: 547.62

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Palbociclib-Succinic acid is a Palbociclib- derivative with a succinic acid linker. The carboxy group of Palbociclib-Succinic acid can be used to conjugate with other molecules such as peptides, proteins, or polymers. Palbociclib-Succinic acid is a useful agent to make Palbociclib-conjugate for drug delivery, nanodrug research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8261 mL	9.1304 mL	18.2608 mL
5 mM	0.3652 mL	1.8261 mL	3.6522 mL
10 mM	0.1826 mL	0.913 mL	1.8261 mL
50 mM	0.0365 mL	0.1826 mL	0.3652 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481