

A-674563

Chemical Properties

CAS No. : 552325-73-2


Formula: C₂₂H₂₂N₄O

Molecular Weight: 358.44

Store at low temperature

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	A-674563 is an Akt1 inhibitor with K_i of 11 nM, modest potent to PKA and >30-fold selective for Akt1 over PKC.
Targets(IC ₅₀)	ERK,Akt,CDK,GSK-3,PKA
In vitro	A-674563 is achieved from A-443654 by replacing the indole with a phenyl moiety and getting oral activity. A-674563 slows proliferation of tumor cells with EC ₅₀ of 0.4 μ M. [1] A-674563 does not inhibit Akt phosphorylation per se, but blocks the phosphorylation of Akt downstream targets in a dose-dependent manner. A-674563 induced Akt blockade results in decreased STS cell downstream target phosphorylation and tumor cell growth inhibition. A-674563 induces G2 cell cycle arrest and apoptosis in STS cells. [2]
In vivo	20 mg/kg A-674563 increases plasma insulin in an oral glucose tolerance test. A-674563 shows no significant monotherapy tumor inhibitory activity; the efficacy of the combination therapy is significantly improved compared to paclitaxel monotherapy. [1] A674563-treated (20 mg/kg/bid, p.o.) mice exhibits slower tumor growth and more than 50% decrease in the tumor volume at the termination of the study compared with that in control group. [2] A-674563 is identified to have drastically improved PK profile with oral bioavailability of 67% in mouse, but is 70-fold less active than A-443654. [3]
Kinase Assay	Akt Kinase Assay: The kinase assay uses His-Akt1 and a biotinylated mouse Bad peptide as substrate. The kinase assay is carried out at room temperature for 30 minutes in 50 μ L of reaction buffer [20 mM HEPES, pH 7.5, 10 mM MgCl ₂ , 0.1% (w/v) Triton X-100, 5 μ M ATP (K_m = 40 μ M), 5 μ M peptide (K_m = 15 μ M), 1 mM DTT, 60 ng of Akt1, and 0.5 μ Ci of [γ - ³³ P]ATP] in the presence of different concentrations of A-674563. Each reaction is stopped by adding 50 μ L of termination buffer (0.1 M EDTA, pH 8.0, and 4 M NaCl). The biotinylated Bad peptides are immobilized on streptavidin-coated FLASH plates. After being washed with PBS-Tween 20 (0.05%), the ³³ P phosphopeptide captured on the FLASH plates is measured with a TopCount Packard Instruments γ counter.
Cell Research	The cells on 96-well plates are gently washed with 200 μ L of PBS. Alamar Blue reagent is diluted 1:10 in normal growth media. The diluted Alamar Blue reagent (100 M) is added to each well on the 96-well plates and incubated until the reaction is complete as per manufacturer's instructions. Analysis is done using an fmax Fluorescence Microplate Reader, set at the excitation wavelength of 544 nm and emission wavelength of 595 nm.

A DRUG SCREENING EXPERT

Cell Research	Data are analyzed using SOFTmax PRO software provided by the manufacturer. (Only for Reference)
---------------	---

Solubility Information

Solubility	H2O: 66 mg/mL (184.13 mM),Sonication is recommended. DMSO: 67 mg/mL (186.92 mM),Sonication is recommended. Ethanol: 15 mg/mL (41.85 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (5.58 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7899 mL	13.9493 mL	27.8987 mL
5 mM	0.558 mL	2.7899 mL	5.5797 mL
10 mM	0.279 mL	1.3949 mL	2.7899 mL
50 mM	0.0558 mL	0.279 mL	0.558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Luo Y, et al, Mol Cancer Ther, 2005, 4(6), 977-986.
Zhu QS, et al, Cancer Res, 2008, 68(8), 2895-2903.
Tatsuya Okuzumi, et al, Mol Biosyst, 2010, 6(8), 1389-1402.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481