


2-[(4-Phenoxyphenyl)formamido]acetic acid

Chemical Properties

CAS No. :	157396-29-7
Formula:	C15H13NO4
Molecular Weight:	271.268
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-[(4-Phenoxyphenyl)formamido]acetic acid, with CAS No. 157396-29-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-[(4-Phenoxyphenyl)formamido]acetic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6864 mL	18.4318 mL	36.8636 mL
5 mM	0.7373 mL	3.6864 mL	7.3727 mL
10 mM	0.3686 mL	1.8432 mL	3.6864 mL
50 mM	0.0737 mL	0.3686 mL	0.7373 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481