

N-(Azido-PEG2)-N-Boc-PEG4-Boc

Chemical Properties

CAS No. : 2093153-83-2


Formula: C₂₆H₅₀N₄O₁₀

Molecular Weight: 578.7

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N-(Azido-PEG2)-N-Boc-PEG4-Boc is a PEG-based bifunctional linker used in the synthesis of PROTACs (proteolysis-targeting chimeras).
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs comprise two distinct ligands linked together: one ligand targets an E3 ubiquitin ligase, while the other targets the specific protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade these target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.728 mL	8.6401 mL	17.2801 mL
5 mM	0.3456 mL	1.728 mL	3.456 mL
10 mM	0.1728 mL	0.864 mL	1.728 mL
50 mM	0.0346 mL	0.1728 mL	0.3456 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481