


Tolbutamide Sodium

Chemical Properties

CAS No. : 473-41-6
 Formula: C₁₂H₁₈N₂NaO₃S
 Molecular Weight: 293.34
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	potassium channel blocker
Targets(IC50)	Apoptosis,Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.409 mL	17.0451 mL	34.0901 mL
5 mM	0.6818 mL	3.409 mL	6.818 mL
10 mM	0.3409 mL	1.7045 mL	3.409 mL
50 mM	0.0682 mL	0.3409 mL	0.6818 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481