


AG-17724

Chemical Properties

CAS No. : 884033-66-3
 Formula: C₂₁H₁₆FN₃O₃
 Molecular Weight: 377.37
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|-------------------------------------------------------------------------|
| Description | AG-17724 is a potent peptidyl-prolyl isomerase (PPIase) Pin1 inhibitor. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.6499 mL | 13.2496 mL | 26.4992 mL |
| 5 mM | 0.530 mL | 2.6499 mL | 5.2998 mL |
| 10 mM | 0.265 mL | 1.325 mL | 2.6499 mL |
| 50 mM | 0.053 mL | 0.265 mL | 0.530 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481