

Ono 11006

## Chemical Properties


CAS No. : 73509-46-3

Formula: C<sub>22</sub>H<sub>36</sub>O<sub>3</sub>

Molecular Weight: 348.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ono 11006 is a thromboxane agonist that dose dependently reduced the positive inotropic effect induced by field stimulation.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8693 mL	14.3464 mL	28.6928 mL
5 mM	0.5739 mL	2.8693 mL	5.7386 mL
10 mM	0.2869 mL	1.4346 mL	2.8693 mL
50 mM	0.0574 mL	0.2869 mL	0.5739 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481