


3'-Azido-3'-deoxy-N6-methyladenosine

Chemical Properties

CAS No. :	889126-08-3
Formula:	C ₁₁ H ₁₄ N ₈ O ₃
Molecular Weight:	306.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-Azido-3'-deoxy-N6-methyladenosine is a Nucleoside Derivative - Azido-nucleoside, 6-Modified purine nucleoside;N-Methylated nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.265 mL	16.3249 mL	32.6499 mL
5 mM	0.653 mL	3.265 mL	6.530 mL
10 mM	0.3265 mL	1.6325 mL	3.265 mL
50 mM	0.0653 mL	0.3265 mL	0.653 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481